
## Technical information Productnumber : 404V-1000





## **Product description**

Hot-dip galvanized drip tray with grid suitable for one IBC container. The capacity of the collection container is 1000 litres.

Guideline for storage and temporary storage for drip trays with regard to fire safety, occupational safety and environmental safety in accordance with PGS 15 (publication series of hazardous substances):

A storage facility must be constructed in such a way that leaked or spilled hazardous liquid container is collected. To this end, the collection container must have a capacity of at least 110% of the largest package. However, if 10% of the total content of the packages is more than 110% of the content of the largest package, then the collection capacity must be equal to 10% of the total content of the packages together. The collection container must be sufficiently resistant to the stored liquids.

The collection capacity only applies to liquids. Empty, uncleaned packaging does not count. The collection facility may also be realized in storage. The choice of location for the collection facility may depend on the hazardous liquids stored. Kruizinga.nl is not liable for the type of use!

Please note! Any form of modification is not allowed. If this is the case, all regulations will lapse.

## Specifications

Article arrangement: New Version: for 1x 1000 I IBC Walls: walls closed Top: 1000kg/m2 Colour: colorless Length (mm): 1460 Height (mm): 620 Weight (kg) 194 Type number: 4041000

Scan the QR code for more information

Type: receptacle for IBC Material: steel Floor: with galvanized grid Subgroup: retention basin steel Surface treatment: galvanized Width (mm): 1460 Loading capacity (kg): 1900 Collection volume (ltr): 1100 EAN Code (Gtin) 8719424006148






















KRUIZINGA.FR® Stockage et ressources de transport

